

Neurofilament light polypeptide

Catalog No: #AP86291



Package Size: #AP86291-1 50ug #AP86291-2 100ug #AP86291-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Neurofilament light polypeptide |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | NF-L,68 kDa neurofilament protein,Neurofilament triplet L protein,Nf68,Nfl |
| Accession No. | Uniprot ID:P19527 |
| Uniprot | P19527 |
| GeneID | 83613; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nefl

Note: This product is for in vitro research use only