

Poly [ADP-ribose] polymerase 1

Catalog No: #AP86292



Package Size: #AP86292-1 50ug #AP86292-2 100ug #AP86292-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Poly [ADP-ribose] polymerase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PARP-1,ADP-ribosyltransferase diphtheria toxin-like 1,ARTD1,NAD(+) ADP-ribosyltransferase 1,ADPRT 1,Poly[ADP-ribose] synthase 1,ADPRT
Accession No.	Uniprot ID:P18493
Uniprot	P18493
GeneID	286764;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PARP1

Note: This product is for in vitro research use only