

Procollagen C-endopeptidase enhancer 1

Catalog No: #AP86293



Package Size: #AP86293-1 50ug #AP86293-2 100ug #AP86293-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Procollagen C-endopeptidase enhancer 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	P14,Procollagen COOH-terminal proteinase enhancer 1,PCPE-1,Type 1 procollagen C-proteinase enhancer protein,Type I procollagen COOH-terminal proteinase enhancer,Pcpe1
Accession No.	Uniprot ID:Q61398
Uniprot	Q61398
GeneID	18542;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pcolce

Note: This product is for in vitro research use only