

Oncoprotein-induced transcript 3 protein

Catalog No: #AP86295



Package Size: #AP86295-1 50ug #AP86295-2 100ug #AP86295-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oncoprotein-induced transcript 3 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Liver-specific zona pellucida domain-containing protein,Lzp
Accession No.	Uniprot ID:Q8R4V5
Uniprot	Q8R4V5
GeneID	18302;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Oit3

Note: This product is for in vitro research use only