

Protocadherin gamma-C3

Catalog No: #AP86297



Package Size: #AP86297-1 50ug #AP86297-2 100ug #AP86297-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protocadherin gamma-C3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PCDH-gamma-C3,Protocadherin-2,Protocadherin-43,PC-43,PCDH2
Accession No.	Uniprot ID:Q9UN70
Uniprot	Q9UN70
GeneID	5098;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCDHGC3

Note: This product is for in vitro research use only