

Phosphatidylethanolamine-binding protein 1

Catalog No: #AP86302



Package Size: #AP86302-1 50ug #AP86302-2 100ug #AP86302-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylethanolamine-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PEBP-1,HCNPPp,Raf kinase inhibitor protein,RKIP,PBP
Accession No.	Uniprot ID:Q3YIX4
Uniprot	Q3YIX4
GeneID	477501;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PEBP1

Note: This product is for in vitro research use only