

Gamma-secretase subunit PEN-2

Catalog No: #AP86303



Package Size: #AP86303-1 50ug #AP86303-2 100ug #AP86303-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-secretase subunit PEN-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Liver regeneration-related protein LRRGT00140, Presenilin enhancer protein 2
Accession No.	Uniprot ID:Q6QI68
Uniprot	Q6QI68
GeneID	292788;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Psenen

Note: This product is for in vitro research use only