

2'-5'-oligoadenylate synthase 3

Catalog No: #AP86306



Package Size: #AP86306-1 50ug #AP86306-2 100ug #AP86306-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2'-5'-oligoadenylate synthase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	(2-5')oligo(A) synthase 3,2',5'-oligoadenylate synthetase-like 10,oasl10
Accession No.	Uniprot ID:Q8VI93
Uniprot	Q8VI93
GeneID	246727;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Oas3

Note: This product is for in vitro research use only