

Lipoamide acyltransferase component of branched-chain alpha-keto acid dehydrogenase complex, mitochondrial

Catalog No: #AP86307

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86307-1 50ug #AP86307-2 100ug #AP86307-3 1mg

Description

Product Name	Lipoamide acyltransferase component of branched-chain alpha-keto acid dehydrogenase complex, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	52 kDa mitochondrial autoantigen of primary biliary cirrhosis, Branched chain 2-oxo-acid dehydrogenase complex component E2, BCOADC-E2, Branched-chain alpha-keto acid dehydrogenase complex component E2, BCKAD-E2, Dihydrolipoamide acetyltransferase component of branched-chain alpha-keto acid dehydrogenase complex, Dihydrolipoamide branched chain transacylase, Dihydrolipoalysine-residue (2-methylpropanoyl)transferase, BCATE2
Accession No.	Uniprot ID:P11182
Uniprot	P11182
GeneID	1629;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: DBT

Note: This product is for in vitro research use only