

NAD(P)H-hydrate epimerase

Catalog No: #AP86309



Package Size: #AP86309-1 50ug #AP86309-2 100ug #AP86309-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NAD(P)H-hydrate epimerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	AI-BP,AIBP
Accession No.	Uniprot ID:Q0PIT9
Uniprot	Q0PIT9
GeneID	780405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NAXE

Note: This product is for in vitro research use only