

Nucleotide-binding oligomerization domain-containing protein 1



Catalog No: #AP86310

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86310-1 50ug #AP86310-2 100ug #AP86310-3 1mg

Description

Product Name	Nucleotide-binding oligomerization domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Caspase recruitment domain-containing protein 4,Card4
Accession No.	Uniprot ID:Q8BHB0
Uniprot	Q8BHB0
GeneID	107607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nod1

Note: This product is for in vitro research use only