

N-fatty-acyl-amino acid synthase/hydrolase PM20D1

Catalog No: #AP86317



Package Size: #AP86317-1 50ug #AP86317-2 100ug #AP86317-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-fatty-acyl-amino acid synthase/hydrolase PM20D1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Peptidase M20 domain-containing protein 1
Accession No.	Uniprot ID:Q2T9M7
Uniprot	Q2T9M7
GeneID	513852;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PM20D1

Note: This product is for in vitro research use only