

Proteinase-activated receptor 1

Catalog No: #AP86318



Package Size: #AP86318-1 50ug #AP86318-2 100ug #AP86318-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteinase-activated receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PAR-1,Coagulation factor II receptor,Thrombin receptor,CF2R,PAR1,TR
Accession No.	Uniprot ID:P25116
Uniprot	P25116
GeneID	2149;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:F2R

Note: This product is for in vitro research use only