

# Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit

Catalog No: #AP86325

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86325-1 50ug #AP86325-2 100ug #AP86325-3 1mg

## Description

Product Name	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	DDOST 48 kDa subunit,OST 50 kDa subunit,OST48
Accession No.	Uniprot ID:P48440
Uniprot	P48440
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:DDOST

Note: This product is for in vitro research use only