

Neprilysin

Catalog No: #AP86326



Package Size: #AP86326-1 50ug #AP86326-2 100ug #AP86326-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neprilysin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Atriopeptidase,Enkephalinase,Neutral endopeptidase 24.11,NEP,Skin fibroblast elastase,SFE,CD10
Accession No.	Uniprot ID:P07861
Uniprot	P07861
GeneID	24590;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mme

Note: This product is for in vitro research use only