

## Pancreas/duodenum homeobox protein 1

Catalog No: #AP86328



Package Size: #AP86328-1 50ug #AP86328-2 100ug #AP86328-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pancreas/duodenum homeobox protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Insulin promoter factor 1,IPF-1,Islet/duodenum homeobox 1,IDX-1,Somatostatin-transactivating factor 1,STF-1,lpf1
Accession No.	Uniprot ID:P52947
Uniprot	P52947
GeneID	29535;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pdx1

Note: This product is for in vitro research use only