

Proenkephalin-A

Catalog No: #AP86329



Package Size: #AP86329-1 50ug #AP86329-2 100ug #AP86329-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proenkephalin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Penk-rs, Penk1
Accession No.	Uniprot ID:P04094
Uniprot	P04094
GeneID	29237;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Penk

Note: This product is for in vitro research use only