

C-type mannose receptor 2

Catalog No: #AP86331



Package Size: #AP86331-1 50ug #AP86331-2 100ug #AP86331-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-type mannose receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Endocytic receptor 180,Macrophage mannose receptor 2,CD280,Endo180
Accession No.	Uniprot ID:Q4TU93
Uniprot	Q4TU93
GeneID	498011;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mrc2

Note: This product is for in vitro research use only