

Metallothionein-1

Catalog No: #AP86332



Package Size: #AP86332-1 50ug #AP86332-2 100ug #AP86332-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metallothionein-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MT-1,Metallothionein-I,MT-I
Accession No.	Uniprot ID:P02802
Uniprot	P02802
GeneID	17748;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mt1

Note: This product is for in vitro research use only