

Metallothionein-2

Catalog No: #AP86333



Package Size: #AP86333-1 50ug #AP86333-2 100ug #AP86333-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Metallothionein-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | MT-2,Metallothionein-2A,Metallothionein-II,MT-II,CES1,MT2 |
| Accession No. | Uniprot ID:P02795 |
| Uniprot | P02795 |
| GeneID | 4502; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MT2A

Note: This product is for in vitro research use only