

## Matrix metalloproteinase-14

Catalog No: #AP86334



Package Size: #AP86334-1 50ug #AP86334-2 100ug #AP86334-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Matrix metalloproteinase-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MMP-14,MMP-X1,MT-MMP,Membrane-type matrix metalloproteinase 1,MT-MMP 1,Membrane-type-1 matrix metalloproteinase,MT1-MMP,Mtmmp
Accession No.	Uniprot ID:P53690
Uniprot	P53690
GeneID	17387;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mmp14

Note: This product is for in vitro research use only