

Phosphoglucomutase-1

Catalog No: #AP86335



Package Size: #AP86335-1 50ug #AP86335-2 100ug #AP86335-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphoglucomutase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PGM 1,Glucose phosphomutase 1,Phosphoglucomutase-2,Pgm2
Accession No.	Uniprot ID:Q9D0F9
Uniprot	Q9D0F9
GeneID	72157;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pgm1

Note: This product is for in vitro research use only