

PCNA-associated factor

Catalog No: #AP86336



Package Size: #AP86336-1 50ug #AP86336-2 100ug #AP86336-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | PCNA-associated factor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | PCNA-associated factor of 15 kDa,PAF15,PCNA-clamp-associated factor,PAF |
| Accession No. | Uniprot ID:Q5E9B2 |
| Uniprot | Q5E9B2 |
| GeneID | 540737; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCLAF

Note: This product is for in vitro research use only