

Prostaglandin D2 receptor

Catalog No: #AP86337



Package Size: #AP86337-1 50ug #AP86337-2 100ug #AP86337-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin D2 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PGD receptor, Prostanoid DP receptor
Accession No.	Uniprot ID:Q9R261
Uniprot	Q9R261
GeneID	498475;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ptgdr

Note: This product is for in vitro research use only