

Epididymal secretory protein E1

Catalog No: #AP86338



Package Size: #AP86338-1 50ug #AP86338-2 100ug #AP86338-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Epididymal secretory protein E1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Human epididymis-specific protein 1,He1,Niemann-Pick disease type C2 protein,HE1
Accession No.	Uniprot ID:P61916
Uniprot	P61916
GeneID	10577;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name: NPC2

Note: This product is for in vitro research use only