

NAD(P)H dehydrogenase [quinone] 1

Catalog No: #AP86339



Package Size: #AP86339-1 50ug #AP86339-2 100ug #AP86339-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD(P)H dehydrogenase [quinone] 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Azoreductase,DT-diaphorase,DTD,Menadione reductase,NAD(P)H:quinone oxidoreductase 1,Phylloquinone reductase,Quinone reductase 1,QR1,Dia4,Nmo1,Nmor1
Accession No.	Uniprot ID:Q64669
Uniprot	Q64669
GeneID	18104;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nqo1

Note: This product is for in vitro research use only