

## Pleckstrin homology-like domain family A member 2

Catalog No: #AP86340



Package Size: #AP86340-1 50ug #AP86340-2 100ug #AP86340-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pleckstrin homology-like domain family A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Imprinted in placenta and liver protein,IPL
Accession No.	Uniprot ID:Q0VC85
Uniprot	Q0VC85
GeneID	618810;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PHLDA2

Note: This product is for in vitro research use only