

Platelet-activating factor acetylhydrolase IB subunit gamma



Catalog No: #AP86341

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86341-1 50ug #AP86341-2 100ug #AP86341-3 1mg

Description

Product Name	Platelet-activating factor acetylhydrolase IB subunit gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PAF acetylhydrolase 29 kDa subunit,PAF-AH 29 kDa subunit,PAF-AH subunit gamma,PAFAH subunit gamma,Platelet-activating factor acetylhydrolase alpha 1 subunit,PAF-AH alpha 1,Pafahg
Accession No.	Uniprot ID:O35263
Uniprot	O35263
GeneID	114113;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pafah1b3

Note: This product is for in vitro research use only