

Myelin P2 protein

Catalog No: #AP86344



Package Size: #AP86344-1 50ug #AP86344-2 100ug #AP86344-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myelin P2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Peripheral myelin protein 2
Accession No.	Uniprot ID:P02689
Uniprot	P02689
GeneID	5375;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PMP2

Note: This product is for in vitro research use only