

Cytosolic phospholipase A2

Catalog No: #AP86346



Package Size: #AP86346-1 50ug #AP86346-2 100ug #AP86346-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosolic phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	cPLA2, Phospholipase A2 group IVA, CPLA2, PLA2G4
Accession No.	Uniprot ID:P49147
Uniprot	P49147
GeneID	396394;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLA2G4A

Note: This product is for in vitro research use only