

Ribosyldihyronicotinamide dehydrogenase [quinone]

Catalog No: #AP86347



Package Size: #AP86347-1 50ug #AP86347-2 100ug #AP86347-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosyldihyronicotinamide dehydrogenase [quinone]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NRH dehydrogenase [quinone] 2, NRH:quinone oxidoreductase 2, Quinone reductase 2, QR2, Nmor2
Accession No.	Uniprot ID:Q9JI75
Uniprot	Q9JI75
GeneID	18105;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nqo2

Note: This product is for in vitro research use only