

## Bile acid receptor

Catalog No: #AP86348



Package Size: #AP86348-1 50ug #AP86348-2 100ug #AP86348-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bile acid receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Farnesoid X-activated receptor,Farnesol receptor HRR-1,Nuclear receptor subfamily 1 group H member 4,Retinoid X receptor-interacting protein 14,RXR-interacting protein 14,BAR,FXR,HRR1,RIP14
Accession No.	Uniprot ID:Q96R11
Uniprot	Q96R11
GeneID	9971;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NR1H4

Note: This product is for in vitro research use only