

Nucleobindin-2

Catalog No: #AP86349



Package Size: #AP86349-1 50ug #AP86349-2 100ug #AP86349-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nucleobindin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DNA-binding protein NEFA,Prepronesfatin,Nefa
Accession No.	Uniprot ID:Q9JI85
Uniprot	Q9JI85
GeneID	59295;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nucb2

Note: This product is for in vitro research use only