

PDZ and LIM domain protein 1

Catalog No: #AP86356



Package Size: #AP86356-1 50ug #AP86356-2 100ug #AP86356-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDZ and LIM domain protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-terminal LIM domain protein 1,Elfin,LIM domain protein CLP-36,Clim1
Accession No.	Uniprot ID:O70400
Uniprot	O70400
GeneID	54132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdim1

Note: This product is for in vitro research use only