

Pyridoxal kinase

Catalog No: #AP86357



Package Size: #AP86357-1 50ug #AP86357-2 100ug #AP86357-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pyridoxal kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pyridoxine kinase,C21orf124,C21orf97,PKH,PNK,PRED79
Accession No.	Uniprot ID:O00764
Uniprot	O00764
GeneID	8566;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDXK

Note: This product is for in vitro research use only