

Proteolipid protein 2

Catalog No: #AP86358



Package Size: #AP86358-1 50ug #AP86358-2 100ug #AP86358-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proteolipid protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q6Y1E2
Uniprot	Q6Y1E2
GeneID	399683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLP2

Note: This product is for in vitro research use only