

Lysophosphatidylcholine acyltransferase 2

Catalog No: #AP86359



Package Size: #AP86359-1 50ug #AP86359-2 100ug #AP86359-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysophosphatidylcholine acyltransferase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LPC acyltransferase 2,1-acylglycerol-3-phosphate O-acyltransferase 11,1-AGP acyltransferase 11,1-acylglycerophosphocholine O-acyltransferase,1-alkylglycerophosphocholine O-acetyltransferase,Acetyl-CoA:lyso-platelet-activating factor acetyltransferase,Acetyl-CoA:lyso-PAF acetyltransferase,Acyltransferase-like 1,Aytl1,Aytl1a,Lpcat2a
Accession No.	Uniprot ID:Q8BYI6
Uniprot	Q8BYI6
GeneID	270084;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lpcat2

Note: This product is for in vitro research use only