

Orexin receptor type 2

Catalog No: #AP86362



Package Size: #AP86362-1 50ug #AP86362-2 100ug #AP86362-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orexin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Ox-2-R,Hypocretin receptor type 2
Accession No.	Uniprot ID:O62809
Uniprot	O62809
GeneID	397428;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HCRTR2

Note: This product is for in vitro research use only