

## Prolyl 4-hydroxylase subunit alpha-3

Catalog No: #AP86364



Package Size: #AP86364-1 50ug #AP86364-2 100ug #AP86364-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Prolyl 4-hydroxylase subunit alpha-3   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | 4-PH alpha-3,Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-3,UNQ711/PRO1374 |
| Accession No.      | Uniprot ID:Q7Z4N8  |
| Uniprot            | Q7Z4N8   |
| GeneID             | 283208;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                       |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:P4HA3

Note: This product is for in vitro research use only