

Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial

Catalog No: #AP86366

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86366-1 50ug #AP86366-2 100ug #AP86366-3 1mg

Description

Product Name	Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PDHE1-A type I
Accession No.	Uniprot ID:A7MB35
Uniprot	A7MB35
GeneID	407109;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PDHA1

Note: This product is for in vitro research use only