

## Plasminogen activator inhibitor 2 type A

Catalog No: #AP86367



Package Size: #AP86367-1 50ug #AP86367-2 100ug #AP86367-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Plasminogen activator inhibitor 2 type A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PAI-2,Serpin B2,Pai2,Planh2
Accession No.	Uniprot ID:P29524
Uniprot	P29524
GeneID	60325;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Serpinb2

Note: This product is for in vitro research use only