

## Proteinase-activated receptor 2

Catalog No: #AP86370



Package Size: #AP86370-1 50ug #AP86370-2 100ug #AP86370-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Proteinase-activated receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PAR-2,Coagulation factor II receptor-like 1,Thrombin receptor-like 1,PAR2
Accession No.	Uniprot ID:Q2HJA4
Uniprot	Q2HJA4
GeneID	526525;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:F2RL1

Note: This product is for in vitro research use only