

Plasma cell-induced resident endoplasmic reticulum protein

Catalog No: #AP86372

Package Size: #AP86372-1 50ug #AP86372-2 100ug #AP86372-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasma cell-induced resident endoplasmic reticulum protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Plasma cell-induced resident ER protein,pERp1,Proapoptotic caspase adapter protein,HSPC190
Accession No.	Uniprot ID:Q8WU39
Uniprot	Q8WU39
GeneID	51237;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PACAP

Note: This product is for in vitro research use only