

Prostaglandin G/H synthase 2

Catalog No: #AP86373



Package Size: #AP86373-1 50ug #AP86373-2 100ug #AP86373-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin G/H synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cyclooxygenase-2, COX-2, Glucocorticoid-regulated inflammatory cyclooxygenase, Gripghs, Macrophage activation-associated marker protein P71/73, PES-2, PHS II, Prostaglandin H2 synthase 2, PGH synthase 2, Prostaglandin-endoperoxide synthase 2, TIS10 protein, Cox-2, Cox2, Pghs-b, Tis10
Accession No.	Uniprot ID:Q05769
Uniprot	Q05769
GeneID	19225;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ptgs2

Note: This product is for in vitro research use only