

Nephrin

Catalog No: #AP86375



Package Size: #AP86375-1 50ug #AP86375-2 100ug #AP86375-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nephrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Renal glomerulus-specific cell adhesion receptor,Nphn
Accession No.	Uniprot ID:Q9R044
Uniprot	Q9R044
GeneID	64563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nphs1

Note: This product is for in vitro research use only