

cAMP and cAMP-inhibited cGMP 3',5'-cyclic phosphodiesterase 10A

Catalog No: #AP86378

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86378-1 50ug #AP86378-2 100ug #AP86378-3 1mg

Description

Product Name	cAMP and cAMP-inhibited cGMP 3',5'-cyclic phosphodiesterase 10A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q9QYJ6
Uniprot	Q9QYJ6
GeneID	63885;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pde10a

Note: This product is for in vitro research use only