

Platelet-derived growth factor subunit B

Catalog No: #AP86379



Package Size: #AP86379-1 50ug #AP86379-2 100ug #AP86379-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor subunit B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PDGF subunit B,PDGF-2,Platelet-derived growth factor B chain,Platelet-derived growth factor beta polypeptide,Proto-oncogene c-Sis,Sis
Accession No.	Uniprot ID:P31240
Uniprot	P31240
GeneID	18591;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdgfb

Note: This product is for in vitro research use only