

Poly [ADP-ribose] polymerase 1

Catalog No: #AP86380



Package Size: #AP86380-1 50ug #AP86380-2 100ug #AP86380-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Poly [ADP-ribose] polymerase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PARP-1,ADP-ribosyltransferase diphtheria toxin-like 1,ARTD1,NAD(+) ADP-ribosyltransferase 1,ADPRT 1,Poly[ADP-ribose] synthase 1,msPARP,Adprp,Adprt,Adprt1
Accession No.	Uniprot ID:P11103
Uniprot	P11103
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Parp1

Note: This product is for in vitro research use only