

Prostate androgen-regulated mucin-like protein 1 homolog

Catalog No: #AP86381

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86381-1 50ug #AP86381-2 100ug #AP86381-3 1mg

Description

Product Name	Prostate androgen-regulated mucin-like protein 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PARM-1
Accession No.	Uniprot ID:Q923D3
Uniprot	Q923D3
GeneID	231440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Parm1

Note: This product is for in vitro research use only