

Proliferating cell nuclear antigen

Catalog No: #AP86382



Package Size: #AP86382-1 50ug #AP86382-2 100ug #AP86382-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proliferating cell nuclear antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PCNA,Cyclin
Accession No.	Uniprot ID:P04961
Uniprot	P04961
GeneID	25737;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pcna

Note: This product is for in vitro research use only